


TrackRAD Project Overview

June 25, 2025

MRI-guided linear accelerator (MRI-linac)

- Combine real-time MR imaging with a linear accelerator for radiation treatment
 - Use 0.35T or 1.5T field strength (less than diagnostic MR but good enough for soft tissue contrast)
 - Linac = high energy X ray beams
 - Beam is often perpendicular to magnetic field
- Uses
 - Real-time tumor tracking to account for motion during breathing
 - Adjust X ray beam in real-time based on MR imaging



MR-LINAC TIMELINE (LAGENDIJK, RADONC 2022)

1999

MRI guided radiotherapy: a
MRI based linear accelerator
Legendijk JTW and Bakker CJG
2000 Proc. ESTRO Istanbul
19th Annual Meeting
Rad Onc vol 56, suppl. 1, pp
S1-S25S

invention

2004



design

2005



Initial exp.

2009



1st prototype

2012



2nd prototype

2014



3rd prototype


2018





Clinical



Frame 0, Patient ID: A_001, Site: abdomen, B = 0.35 [T], Fs = 8 [Hz]





Frame: 3, Delta = 2.2 [mm], $(V, V_0) = (14.9, 15.4) \text{ [cm}^2\text{]}$



MR object tracking status quo

- Current clinically available solutions use:
 - Deformable image registration – vector field maps points from frame to frame
 - Template matching – slide a reference image across the target, calculate similarity scores, optimize
- Both struggle with large, non-rigid motion
- Best practice = turn off the beam when “large motion” is detected (gating)
 - This wastes valuable scanner time
- MRI adds further complexity
 - Lack of standardization in field strengths, imaging protocols
 - Variance in image quality and artifacts
 - Pixel intensities are relative and vary b/w scans (vs. CT where 0 intensity == water)




TrackRAD2025

- Goal = improve real-time tumor tracking algorithm for MRI-linac
- Task
 - Input
 - **Sequence of 2D MRI frames (AKA a video), 100-200 frames (~10-30s at 8 fps)**
 - **First frame labeled (all pixels pertaining to the tumor)**
 - Metadata (field strength, frame rate)
 - Output = labels for all subsequent frames
- Metrics = image similarity (Dice coefficient), runtime (<1s/frame)
- Dataset
 - From 6 institutions (3 Dutch, 1 German, 1 Australian, 1 Chinese)
 - 0.35T and 1.5T MRI
 - Labeled = 108 patient cases = 10_000 frames
 - 50 for training, 8 + 50 held out for validation and testing
 - **Unlabeled = 477 patient cases = 2_800_000 frames**


Segment Anything Model (SAM) (2023)

- Foundation vision model by Meta
- Transformer architecture
- Trained on 11M images + 1.1B segmentation masks
- “Promptable” = use positive/negative points, bounding box
- Good “zero-shot” performance = applied to a task without fine-tuning
- Can be fine-tuned (trained with domain-specific data)
- <https://segment-anything.com/demo>



SAM2 (2024)

- Adds video object segmentation (VOS) functionality
 - **Given one frame + label, infer labels in all subsequent frames**
 - Robust to deformation, object occlusion, entering and leaving the field of view
- Trained on video
 - Training data includes the previous + following 8 frames
- <https://ai.meta.com/sam2/>



SAM2 evaluation – it's pretty good!

Zero-shot with 50 labeled test cases

"aggregates": {

"dice_similarity_coefficient": 0.9071, # 2 * image overlap / total pixels, [0,1]

"hausdorff_distance_95": 3.546, # for each point, min dist to any other point, 95th pct

"surface_distance_average": 1.303, # for each point, min dist to any point, mean

"center_distance": 1.595, # dist b/w centers of mass

"relative_d98_dose": 0.9434, # est accuracy of X-ray dose, scatter and stuff


"time_per_frame": 0.1900, # must be <1

}

Modeling plan

1. Use SAM2 for tumor tracking
 - a. Frame sequence + first frame label → Infer subsequent frames
2. Use SAM2 to expedite annotation of unlabeled frame sequences
 - a. Assist in annotating first frames w/bounding box → generate mask
 - b. Infer masks for subsequent frames
 - c. Validate mask quality using traditional rule-based methods
 - d. Flag bad frames for relabeling
3. Fine-tune SAM2 using the labeled data
4. (Maybe) Find-tune with other open-access datasets
5. Deploy fine-tuned model

Data provided by the challenge



[Link](#)

Labeling tool

<https://huggingface.co/spaces/mayo-radonc/bouncing-target>

- HuggingFace = tech company that provides software + services for ML apps
 - Hosts datasets (including TrackRAD)
 - Hosts applications, managing all the hardware
- Unlabeled data = 477 patients, ~2000 sequences, 2.8M frames
- What do I label?
 - Any and all objects
 - Goal = not “find the tumor” but “get better at black and white, fuzzy, medical images”

The importance of labeled data

- Labeling is a critical part of this project
- Anyone can download a foundation model from the internet and run it
- **Our competitive edge = rapidly labeling a large dataset for finetuning**
- 2.8M frames is a lot of data, but our plan is to can automate 99% of it, manually labeling a few key frames per video
- We will also look into other open-access imaging datasets

Up next

1. Complete labeling tool (ETA early July)
2. **Start labeling!**
3. Develop annotation pipeline
 - a. Download labeled frames from HuggingFace
 - b. Use SAM2 to infer labels for subsequent frames
 - c. **Validate frames, flag bad frames**
4. Fine-tune SAM2 with labeled frames in batches
5. **(Maybe) Fine-tune SAM2 with other open source datasets**
6. **Optimize SAM2 parameters (confidence threshold)**
7. **Submit model (August 15)**

What's in it for me?

- First-hand experience with AI tools in medical imaging + computer vision
- The top 5 teams (team limit = 5 people) get
 - A cash prize (\$1000 to \$200, will distribute amongst everyone who participates)
 - Invited to present at MICCAI 2025, 9/23-27 in S. Korea
 - Included in a “challenge” paper submitted to Medical Image Analysis or a similar journal
 - <https://trackrad2025.grand-challenge.org/teams/>
- We're currently in 1st place among preliminary submissions (as of June 25), so it's definitely doable!
 - <https://trackrad2025.grand-challenge.org/evaluation/preliminary-testing/leaderboard/>
- Additional paper(s) as appropriate